

## 10x PBS Solution (Phosphate Buffered Saline)

### Product Description

Phosphate Buffered Saline (PBS) is a balanced salt solution that can be used in a variety of biological assays.

Note: The salt concentration of this buffer is at 10x. A single (1x) strength solution contains 0.137M NaCl, 0.0027M KCl, 0.008M Na<sub>2</sub>HPO<sub>4</sub> and 0.002M KH<sub>2</sub>PO<sub>4</sub>.

**Product No.:** PD8117

**Grade:** Biotech

**Available Volume:** 4L

**Storage Condition:** 18 – 25 °C



BIO BASIC INC.