

SDS, 20% (w/v) Solution

Product Description

SDS is a protein denaturant used in polyacrylamide gel electrophoresis for molecular weight determinations. It is also used to dissociate nucleic acid-protein complexes in DNA extraction protocols, to disrupt cell membranes, and to prepare pre-hybridization and hybridization solutions.

Note: This contains 20% sodium dodecyl sulfate (SDS).

Product No.: SD8119

Grade: Biotech

Available Volume: 100ml, 500ml, 1L

Storage Condition: 18 – 25 °C



BIO BASIC INC.