

TT-SDS (Tris-Tricine-SDS Buffer) Premix Powder

Product Description

Tris-Tricine (TT)-SDS Buffer replaces Glycine with Tricine for use in the standard TG and TG-SDS buffer system. It allows better resolution of low molecular weight proteins and peptides (less than 10 kDa).

Note: A pack makes 1L of 10x TT-SDS Buffer or 10L of 1x TT-SDS Buffer. A single strength (1x) buffer contains 0.01M Tris base, 0.01M Tricine and 0.01% SDS.

Product No.: TD8135

Grade: Biotech

Available Volume: 10L

Storage Condition: 18 – 25 °C

