

TBS Buffer (Tris-Buffered Saline), 20x Solution

Product Description

This product is ideal for use as a wash buffer and antibody diluent for ELISA, Western blotting and other immunoassays.

Note: The salt concentration in this product in this buffer is at 20x. A single (1x) strength is 50 mM Tris-Cl, 150 mM NaCl pH 7.5.

Product No.: UA0045

Grade: Ultra-Pure, 0.22 μ M filtered

Available Volume: 500ml, 1L

Storage Condition: 18 – 25 °C.

