

## TBST (TBS-Tween 20), 10X

### Product Description

This product is commonly used as wash buffer in Western Blot, ELISA, and immunoassays. It can also be used as blocking buffer for plate based assays.

Note: The salt concentration of this buffer is at 10x. A single (1x) strength contains 0.025M Tris, 0.15M NaCl, 0.05% Tween-20, pH 7.5.

**Product No.:** UA0046

**Grade:** Ultra-Pure, 0.22  $\mu$ M filtered

**Available Volume:** 500ml, 1L

**Storage Condition:** 18 – 25 °C. Protect from moisture.



BIO BASIC INC.